



SONDERVERSORGUNG



BERLIN

STOMOCUR® Drainage Protect Bed Bag - DPB2000

Safety and comfort for a restful sleep at night

The high-capacity STOMOCUR® Drainage Protect Bed Bag

Your advantages:

Integrated adapter for connection to all FORLIFE Drainage Protect pouches

130 cm long, large-lumen and kink-resistant tube

Large outlet, also for undigested residues

Integrated suspension

STOMOCUR® Drainage Protect Bed Bag

Pack contains 5 pieces

Tube length	Capacity	Article number
130 cm	2.000 ml	DPB2000

Compatible with:


- STOMOCUR® Drainage Protect
- STOMOCUR® Drainage Protect konvex
- STOMOCUR® Drainage Protect medium konvex
- STOMOCUR® Drainage Protect Clic

Packs contain 10 pieces. CDP70T contains 30 pieces.

Stoma opening	Plateau size	Article number
20 – 60 mm	-	DP5020
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3
20 – 47 mm	55 mm	DP20M1
15 – 40 mm	45 mm	DP15M2
10 – 28 mm	35 mm	DP10M3
Locking ring: 70 mm	-	CDP70T



1 The integrated adapter allows it to be easily connected to all Drainage Protect pouches from FORLIFE.



2 The 130 cm long kink-resistant tube is passable for larger, undigested residues.



3 The outlet with its large diameter allows easy emptying and cleaning.



4 The integrated suspension device makes the bag easy to place on the bed.

