

PRIMO® 

Reliability has a new colour

Discover our new range
of STOMOCUR® Primo®
pouches. Discreet protection
and noticeable comfort.


by
FORLIFE
BERLIN

Empowering yourself.

Our trusted skin protection with an ultra soft black pouch.

Our new Primo® product range is a combination of the proven benefits of our successful STOMOCUR® Convex and Medium Convex series with a modern, stylish pouch design and a high-quality textile cover.







Alginat

With its natural, supportive abilities to improve and protect skin health through optimal moisture absorption, it offers preventative advantages.

Combiflange Hydrocolloid

Our skin-friendly Combiflange Hydrocolloid skin protection material is soft, flexible and secure to maintain healthy skin, while its tapered edge helps reduce leaks.





Ultimate softness

Soft, smooth textile cover



Blackout weave

100% Opaque



Hydrophobic surface

Water-repellent









Designed for comfort

Asymmetrical shape



Available convexities within the Primo® range

Our convex and medium convex baseplates are contoured and specially developed to adapt flexibly to the contours of your body and thus ensure optimum wearing comfort.





Ileostomy pouch size 155 x 292 mm, capacity 460 ml, 10 pcs.

Stoma Opening	Plateau Size	Article Number
15 - 43 mm	50 mm	IAH15V1B
15 - 35 mm	40 mm	IAH15V2B
10 - 25 mm	30 mm	IAH10V3B

Stoma Opening	Plateau Size	Article Number
15 - 43 mm	50 mm	IAA15V1B
15 - 35 mm	40 mm	IAA15V2B
10 - 25 mm	30 mm	IAA10V3B

Colostomy pouch size 154 x 202 mm, capacity 320 ml, 30 pcs.

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	CAH20V1B
15 - 35 mm	40 mm	CAH15V2B
10 - 25 mm	30 mm	CAH10V3B

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	CAA20V1B
15 - 35 mm	40 mm	CAA15V2B
10 - 25 mm	30 mm	CAA10V3B



Ileostomy pouch size 155 x 292 mm, capacity 460 ml, 10 pcs.

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IAH15M1B
15 - 40 mm	45 mm	IAH15M2B
10 - 28 mm	35 mm	IAH10M3B

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IAA15M1B
15 - 40 mm	45 mm	IAA15M2B
10 - 28 mm	35 mm	IAA10M3B

Colostomy pouch size 154 x 202 mm, capacity 320 ml, 30 pcs.

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CAH20M1B
15 - 40 mm	45 mm	CAH15M2B
10 - 28 mm	35 mm	CAH10M3B

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CAA20M1B
15 - 40 mm	45 mm	CAA15M2B
10 - 28 mm	35 mm	CAA10M3B



For Life
Produktions- und Vertriebsgesellschaft
für Heil- und Hilfsmittel mbH

Am Studio 16
12489 Berlin, Germany

mail@forlife.info
www.forlife.info



FORLIFE
BERLIN

Empowering yourself.