


We are expanding our range of ostomy caps.
STOMOCUR® Select Cap now in three sizes with AND without filter!

SAFE, DISCREET AND RELIABLE.



Your benefits:

Wide range of different ostomy caps

Hypoallergenic Hydrocolloid skin protection

Very good adhesion and skin compatibility

Absorbent insert for more safety

With soft non-woven on both sides

Easy to use



STOMOCUR® Select Cap

Package contents 30 pieces

	Stoma opening	Article number	Cap size	Adhesive surface Ø	Filter
NEW	20 - 40 mm	CSC20SNF	60 mm	50 mm	No
	20 - 50 mm	CSC20	75 mm	60 mm	Yes
NEW	25 - 50 mm	CSC25NF	75 mm	60 mm	No
NEW	30 - 60 mm	CSC30LNF	85 mm	70 mm	No

What does the article number tell us?

Position 1: Pouch	Position 2: Opening	Position 3: Supplement 1	Position 4: Supplement 2
CSC = Colostomy Select Cap	20, 25, 30	L = Large S = Small	NF = No Filter

