



**STOMOCUR<sup>®</sup>**

Product Guide 2025  
Stomacare

## Introduction

**FORLIFE has stood for high-quality products for over 30 years, enabling economic activities for medical retailers, clinics and pharmacies while maintaining the highest quality standards – Made in Berlin.**

We have deliberately chosen to conduct research, development and production exclusively in Germany. As well as trusting in the extensive experience of our employees, the location is also part of our corporate philosophy. We are delighted that we have successfully brought together long-standing FORLIFE employees with younger talents, combining proven technologies with innovation - ensuring long-term economic production of high-quality products.

Our goal is to provide ideal care for the stoma, allowing ostomates to get on with their lives in comfort and security. This is the objective of our products and our daily work.

Your FORLIFE Team

For you. For me. For Life

*Empowering yourself.*

We help people to live a self-empowered and carefree life.





## Primo



Our pouches with ultra-soft, 100% opaque and water-repellent cover in an elegant black design and asymmetrical pouch shape.

Page 8



## Hydrocolloid

Our hydrocolloid supply offers superior properties – soft, flexible and reliable.

Page 13



## Combiflange Hydrocolloid

Our Combiflange Hydrocolloid supply offers superior properties with an outstanding extra. Two hydrocolloid layers bonded together ensuring even better adhesion and safety.

Page 13



## Convex

7 mm convexity.

Our convex products adapt specifically to the contours of your body and thus ensure optimum wearing comfort.

Page 18



## Medium Convex

3,5 mm convexity.

These products exert a medium amount of pressure on the skin around the stoma preventing leaks and skin irritations as well as increasing comfort.

Page 22



## Alginate

An alginate-containing skin protection. Promotes skin integrity, absorbs moisture, and aids in healing damaged skin.

Page 25



## Clic System

A versatile, flexible and discreet two-piece system.

Page 28



## Urostomy

Our extensive urostomy range is well designed. With two different outlet options, personal preferences can be met flexible.

Page 33



## Microskin

Our thinnest and most flexible skin protection.

Page 40



## Gelatine Free

Using a hydrocolloid skin protection without gelatine or other animal components. Gentle on the skin.

Page 43



## Special Supply

Wide range for very unique medical requirements.

Page 45



## Accessories

A comprehensive range of accessories.


Page 51

## Always well supplied with STOMOCUR®


Our STOMOCUR® pouches are developed and manufactured in Berlin using modern production technologies and high-quality materials. They are characterised by innovative and unique properties that meet medical requirements of healthcare professionals and patients.

### Pouch Shapes

Symmetric



Asymmetric



Outlets



Soft outlet -  
Urostomy

"Twist and Safe"  
tap valve -  
Urostomy

"Safe and Easy"  
closure -  
Ileostomy


Large soft outlet -  
Drainage

### Systems

One-piece



Two-piece




### Cover Material


Non-woven cover



Textile cover



NEW



### "Window"-Technology



Many of our pouches are featuring a window invented by FORLIFE. This enables easy positioning and observation.

### Filter




Activated carbon filter – Colostomy

MPC filter – Ileostomy

MPC filter technology (Moisture Protect Composition) = additional protective layer that prevents liquid absorption.

### Convexities



Further details on the respective product ranges can be found on the following pages of our STOMOCUR® product guide.

Primo

Hydrocolloid

Convex

Medium Convex

Alginate

Clic System

Urostomy

Gelatine Free

Special Supply

Accessories

# PRIMO®

NEW

Our new Primo® product range is a combination of the proven benefits of our successful STOMOCUR® Convex and Medium Convex series, with a modern, stylish pouch design and a high-quality textile cover.

Choose between our Alginate and Combiflange Hydrocolloid skin protection and be assured of our safe, skin-friendly and high-performance products – specially developed to adapt flexibly and comfortably to your everyday life.

**Features:**

- High level of comfort
- Easy to use
- Asymmetrical shape
- Elegant design
- 100% opaque

**Benefits:**

- Ultra soft textile cover on both sides
- Window technology
- Discreet design in neutral black
- Water-repellent surface
- With proven "Safe and Easy" closure for ileostomy pouches including black closure components
- Maintains convexity during use and does not flatten against the skin
- Integrated belt attachment




  
 Ultra-soft cover

Available in 3 plateau sizes




Available soon!

**STOMOCUR® Primo® konvex**

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 118 mm (Ø), pouch size 154 x 202 mm, capacity 320 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	CAH20V1B
15 - 35 mm	40 mm	CAH15V2B
10 - 25 mm	30 mm	CAH10V3B



Easy to apply

**STOMOCUR® Primo® konvex**

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 155 x 292 mm, capacity 460 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 43 mm	50 mm	IAH15V1B
15 - 35 mm	40 mm	IAH15V2B
10 - 25 mm	30 mm	IAH10V3B



100% opaque



### STOMOCUR® Primo® Alginat konvex


Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 118 mm (Ø), pouch size 154 x 202 mm, capacity 320 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	CAA20V1B
15 - 35 mm	40 mm	CAA15V2B
10 - 25 mm	30 mm	CAA10V3B

Available soon!

Easy to apply

Water-repellent



### STOMOCUR® Primo® medium konvex

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 118 mm (Ø), pouch size 154 x 202 mm, capacity 320 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CAH20M1B
15 - 40 mm	45 mm	CAH15M2B
10 - 28 mm	35 mm	CAH10M3B

Available soon!



### STOMOCUR® Primo® Alginat konvex

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 155 x 292 mm, capacity 460 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 43 mm	50 mm	IAA15V1B
15 - 35 mm	40 mm	IAA15V2B
10 - 25 mm	30 mm	IAA10V3B

For sensitive skin




### STOMOCUR® Primo® medium konvex

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 155 x 292 mm, capacity 460 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IAH15M1B
15 - 40 mm	45 mm	IAH15M2B
10 - 28 mm	35 mm	IAH10M3B

Easy to apply

Ultra-soft cover



**STOMOCUR® Primo®**  
Alginat medium konvex

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 118 mm (Ø), pouch size 154 x 202 mm, capacity 320 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CAA20M1B
15 - 40 mm	45 mm	CAA15M2B
10 - 28 mm	35 mm	CAA10M3B

Available soon!

For sensitive skin



**STOMOCUR® Primo®**  
Alginat medium konvex

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 155 x 292 mm, capacity 460 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IAA15M1B
15 - 40 mm	45 mm	IAA15M2B
10 - 28 mm	35 mm	IAA10M3B

For sensitive skin



**HYDROKOLLOID**



**COMBIFLANGE**



**How does Hydrocolloid work?**

Hydrocolloid adheres to the skin through the absorption of heat and moisture. This moisture absorption causes the baseplate to swell slightly, ensuring good adhesion to the skin and preventing leaks and skin irritation.

**Features:**

- Hypoallergenic adhesive
- Gentle on the skin
- High level of comfort
- Easy to use
- Highly absorbent adhesive
- High tack adhesion

**Advantages of Combiflange Hydrocolloid in comparison to Standard Hydrocolloid:**

- 0,4 mm thicker than Standard Hydrocolloid
- 10-times higher adhesive strength compared to Standard Hydrocolloid
- Two hydrocolloid layers bonded together ensure even better absorption of the exudating skin, as well as an optimal seal around the stoma
- Available across the range in a variety of one-piece convex products, flat urostomy pouches as well as two-piece convex and flat baseplates
- Thin, tapered border: ensuring a good secure fit even on uneven skin and provides great flexibility



**STOMOCUR® Select large**

Window, filter, soft cover both sides, adhesive surface 135 × 135 mm, pouch size 160 × 220 mm, capacity 470 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 90 mm	CSH20L
25 mm	CSH25L
30 mm	CSH30L
35 mm	CSH35L
40 mm	CSH40L

Cuts up to 90 mm



**STOMOCUR® Select small**

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 130 × 190 mm, capacity 250 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSH20S
25 mm	CSH25S
30 mm	CSH30S
35 mm	CSH35S
40 mm	CSH40S



**STOMOCUR® Select**

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 145 × 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSH20
25 mm	CSH25
30 mm	CSH30
35 mm	CSH35
40 mm	CSH40



**STOMOCUR® Select extra small**

Window, filter, soft cover both sides, adhesive surface 80 mm (Ø), pouch size 100 × 160 mm, capacity 140 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 50 mm	CSH10XS

Also suitable for child care

**STOMOCUR® Select extra small transparent**

Transparent, filter, soft cover one side

Stoma Opening	Article Number
10 – 50 mm	CSH10XST





### STOMOCUR® Protect

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 110 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 70 mm	IPH2070
15 – 60 mm	IPH15
20 mm	IPH20
25 mm	IPH25
30 mm	IPH30
35 mm	IPH35
40 mm	IPH40

### STOMOCUR® Protect transparent

Transparent, filter, soft cover one side, "Safe and Easy" closure

Stoma Opening	Article Number
15 – 60 mm	IPH15T

### STOMOCUR® Protect small

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 110 mm (Ø), pouch size 150 x 190 mm, capacity 400 ml, package contents 30 pieces

Stoma Opening	Article Number
15 – 60 mm	IPH15S



### STOMOCUR® Select Cap

Package contents 30 pieces

Stoma Opening	Article Number	Cap Size	Adhesive Surface Ø	Filter
20 – 40 mm	CSC20SNF	60 mm	50 mm	No
20 – 50 mm	CSC20	75 mm	60 mm	Yes
25 – 50 mm	CSC25NF	75 mm	60 mm	No
30 – 60 mm	CSC30LNF	85 mm	70 mm	No



### STOMOCUR® Drainage Protect

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Article Number
20 – 60 mm	DP5020

Our convex products adapt specifically to the contours of your body and thus ensure optimum wearing comfort. Available as Combiflange Hydrocolloid and Alginate options.



**Features:**

- Hypoallergenic adhesive
- High tack adhesion
- Tapered edge with integral border
- Superior leakage protection
- Soft and flexible
- Anatomically shaped

**Benefits:**

- 7 mm height of convexity
- Four convexities
- Comprehensive range of cut-to-fit and pre-cut sizes
- Efficient filter system for odour reduction
- Integrated belt attachment

**Convex supply:**



**STOMOCUR® Select konvex**

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 57 mm	60 mm	CSH20V0
20 - 43 mm	50 mm	CSH20V1
25 mm		CSH25V1
30 mm		CSH30V1
35 mm		CSH35V1
40 mm		CSH40V1
15 - 35 mm	40 mm	CSH15V2
20 mm		CSH20V2
25 mm		CSH25V2
30 mm		CSH30V2
10 - 25 mm	30 mm	CSH10V3
15 mm		CSH15V3
20 mm		CSH20V3




**STOMOCUR® Protect konvex**

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 57 mm	60 mm	IPH20V0
15 - 43 mm	50 mm	IPH15V1
20 mm		IPH20V1
25 mm		IPH25V1
30 mm		IPH30V1
35 mm		IPH35V1
40 mm	40 mm	IPH40V1
15 - 35 mm		IPH15V2
20 mm		IPH20V2
25 mm		IPH25V2
30 mm	30 mm	IPH30V2
10 - 25 mm		IPH10V3
15 mm		IPH15V3
20 mm		IPH20V3

Easy to apply



For sensitive skin

**STOMOCUR® Select Alginat konvex**

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 × 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 – 43 mm	50 mm	CSA20V1
25 mm		CSA25V1
30 mm		CSA30V1
35 mm		CSA35V1
40 mm		CSA40V1
15 – 35 mm	40 mm	CSA15V2
20 mm		CSA20V2
25 mm		CSA25V2
30 mm		CSA30V2
10 – 25 mm	30 mm	CSA10V3
15 mm		CSA15V3
20 mm		CSA20V3



Soft and flexible

**STOMOCUR® Drainage Protect konvex**

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 × 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3



For sensitive skin

**STOMOCUR® Protect Alginat konvex**

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 150 × 220 mm, capacity 30 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 – 43 mm	50 mm	IPA15V1
20 mm		IPA20V1
25 mm		IPA25V1
30 mm		IPA30V1
35 mm		IPA35V1
40 mm		IPA40V1
15 – 35 mm	40 mm	IPA15V2
10 – 25 mm	30 mm	IPA10V3

Our Medium Konvex range offers the same outstanding flexibility and comfort as our Konvex range, but slightly shallower with a convex height of 3,5 mm. Available as Combiflange Hydrocolloid.



### Features:

- Hypoallergenic adhesive
- High tack adhesion
- Tapered edge with integral border
- Superior leakage protection
- Soft and flexible
- Anatomically shaped

### Benefits:

- 3,5 mm height of convexity
- Three convexities
- Simpler alternative to using a seal with a flat pouch
- Comprehensive range of cut-to-fit and pre-cut baseplates
- Efficient filter system for odour reduction
- Integrated belt attachment

### Medium convex supply:



### STOMOCUR® Select medium konvex

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CSH20M1
25 mm		CSH25M1
30 mm		CSH30M1
35 mm		CSH35M1
40 mm		CSH40M1
15 - 40 mm	45 mm	CSH15M2
20 mm		CSH20M2
25 mm		CSH25M2
10 - 28 mm	35 mm	CSH10M3



### STOMOCUR® Protect medium konvex

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IPH15M1
20 mm		IPH20M1
25 mm		IPH25M1
30 mm		IPH30M1
35 mm		IPH35M1
40 mm	45 mm	IPH40M1
15 - 40 mm		IPH15M2
20 mm		IPH20M2
25 mm	IPH25M2	
10 - 28 mm	35 mm	IPH10M3

### STOMOCUR® Protect medium konvex transparent

Transparent, filter, soft cover one side, "Safe and Easy" closure

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IPH15M1T
15 - 40 mm	45 mm	IPH15M2T
10 - 28 mm	35 mm	IPH10M3T



**STOMOCUR® Drainage Protect**  
medium konvex

Window, filter, soft cover both sides, large soft outlet,  
adhesive surface 118 mm (Ø), pouch size 150 × 220 mm,  
capacity 540 ml, package contents 10 pieces



Soft and  
flexible

Stoma Opening	Plateau Size	Article Number
20 – 47 mm	55 mm	DP20M1
15 – 40 mm	45 mm	DP15M2
10 – 28 mm	35 mm	DP10M3



Alginate is a naturally occurring, bioactive hydrocolloid obtained from brown seaweed. Due to its positive properties in enhancing skin health, alginate has been used to treat affected skin for many years.

Available in a variety of sizes for both drainable and closed pouches.

**How does Alginate work?**

Mixed with water, it creates an elastic gel that accurately adapts to any surface conditions. It can absorb up to 20 times its own weight in liquid and is therefore ideal to enhance wound healing.

**Features:**

- Promotes healing and healthy skin
- Highly absorbent adhesive
- Anatomically shaped
- High level of comfort
- Perfect for skin that is prone to allergies

**Benefits:**

- Contains 15% alginate: This is significant enough to promote skin healing
- Large adhesive area
- Available in two different types of filters: MPC-Filter (Moisture Protect Composition) for ileostomy pouches and standard activated carbon filter for colostomy pouches







For sensitive skin

### STOMOCUR® Select Alginat

Window, filter, soft cover both sides, adhesive surface 110 × 110mm, pouch size 145 × 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSA20
25 mm	CSA25
30 mm	CSA30
35 mm	CSA35
40 mm	CSA40

Easy to apply



For sensitive skin

### STOMOCUR® Protect Alginat

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 110 mm (Ø), pouch size 150 × 220 mm, capacity 530 ml, package contents 30 pieces

Stoma Opening	Article Number
15 – 60 mm	IPA15
20 mm	IPA20
25 mm	IPA25
30 mm	IPA30
35 mm	IPA35
40 mm	IPA40



For sensitive skin

### STOMOCUR® Select Alginat small

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 130 × 190 mm, capacity 250 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSA20S
25 mm	CSA25S
30 mm	CSA30S
35 mm	CSA35S
40 mm	CSA40S

Convex colostomy pouches with Alginate skin protection and black cover can be found on page 10 & 12, with beige cover on page 20.

Convex ileostomy pouches with Alginate skin protection and black cover can be found on page 10 & 12, with beige cover on page 20.



Looking for even more flexibility? Our Two-Piece Clic System makes it easy to change pouches without removing the baseplate. Available in a wide range of flat and convex baseplates for closed, drainable, urostomy and drainage pouches in different locking ring sizes.

Available as Combiflange Hydrocolloid and Alginate skin protection.

**Features:**

- Gentle on the skin
- Locking ring that can gripped behind
- Audible "click"
- Easy to use

**Benefits:**

- 7 mm height of convexity
- Four convexities
- Three different ring sizes (45 mm, 57 mm and 70 mm)
- Comprehensive range of stoma openings, available in cut to fit and pre-cut baseplates
- Integrated belt attachment



Two-Piece Supply | Clic System | Baseplates



**STOMOCUR® Clic**

Plate size 118 mm (Ø), package contents 5 pieces

	Stoma Opening	Article Number
45	15 – 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
57	15 – 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
70	15 – 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050



**STOMOCUR® Clic Alginat**

Plate size 118 mm (Ø), package contents 5 pieces

	Stoma Opening	Article Number
45	15 – 35 mm	CA4515
57	15 – 40 mm	CA5715



Available in 3 plateau sizes

### STOMOCUR® Clic konvex

Plate size 118 mm (Ø), package contents 3 pieces

	Stoma Opening	Plateau Size	Article Number
45	15 – 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
57	15 – 25 mm	30 mm	CK4515V3
	15 – 43 mm	50 mm	CK5715V1
20 mm	CK5720V1		
25 mm	CK5725V1		
30 mm	CK5730V1		
35 mm	CK5735V1		
	15 – 35 mm	40 mm	CK5715V2
	15 – 25 mm	30 mm	CK5715V3



### STOMOCUR® Select Clic

Window, filter, soft cover both sides, pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CC45
57	57 mm	CC57
70	70 mm	CC70



### STOMOCUR® Clic Alginat konvex

Plate size 118 mm (Ø), package contents 3 pieces

	Stoma Opening	Plateau Size	Article Number
45	15 – 35 mm	40 mm	CA4515V2
	15 – 25 mm	30 mm	CA4515V3
	15 – 43 mm	50 mm	CA5715V1
57	15 – 35 mm	40 mm	CA5715V2
	15 – 25 mm	30 mm	CA5715V3



Soft and discreet closure

### STOMOCUR® Protect Clic

Window, filter, soft cover both sides, "Safe and Easy" closure, pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CI45
57	57 mm	CI57
70	70 mm	CI70



Soft and discreet closure

**STOMOCUR® Protect Clic small**

Window, filter, soft cover both sides, "Safe and Easy" closure, pouch size 150 x 190 mm, capacity 400 ml, package contents 30 pieces


Locking Ring Size	Article Number
45	CI45S



**STOMOCUR® Drainage Protect Clic transparent**

Transparent, filter, soft cover one side, large soft outlet, pouch size 150 x 220 mm, capacity 540 ml, package contents 30 pieces

Locking Ring Size	Article Number
70	CDP70T



Our extensive urostomy range is designed with a unique anti-twist shape - makes emptying easier. Available with Combiflange Hydrocolloid, Microskin and gelatine free skin protection in flat, convex or medium convex.

Also for this type of supply, our range offers two-piece solutions.

**Features:**

- Tapered edge with integral border
- Easy to use
- Thin and flexible
- Anatomically shaped
- High level of comfort

**Benefits:**


- Optional 'Safe' soft bung with soft outlet or "Twist and Safe" tap valve
- Unique anti-twist shape
- Drip cap for "Twist and Safe" tap valve
- Non-return valve stops urine back flowing over stoma
- Easy connection to drainage supplies
- Also available with a transparent pouch back



**Microskin - The special skin protection for urostomy supplies**

Due to the particularly thin and flexible skin protection layer, the product adapts perfectly to the skin and thus offers secure adhesion and flexibility for maximum freedom of movement. The Microskin skin protection is permeable to water vapour and therefore prevents unpleasant sweating.





Soft and comfortable

### STOMOCUR® Safe

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 × 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 60 mm	USH10F
15 mm	USH15F
20 mm	USH20F
25 mm	USH25F
30 mm	USH30F

### STOMOCUR® Safe transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Article Number
10 – 60 mm	USH10FT



"Twist and Safe" tap valve

### STOMOCUR® Twist & Safe


Window, soft cover both sides, tap valve, adhesive surface 118 mm (Ø), pouch size 130 × 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 60 mm	USH10H
15 mm	USH15H
20 mm	USH20H
25 mm	USH25H
30 mm	USH30H

### STOMOCUR® Twist & Safe transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Article Number
10 – 60 mm	USH10HT



Soft and comfortable

### STOMOCUR® Safe konvex

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 × 195 mm, capacity 420 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
10 – 43 mm	50 mm	USH10V1F
15 mm		USH15V1F
20 mm		USH20V1F
25 mm		USH25V1F
30 mm		USH30V1F
35 mm		USH35V1F
10 – 35 mm	40 mm	USH10V2F
15 mm		USH15V2F
20 mm		USH20V2F
25 mm		USH25V2F
30 mm		USH30V2F
10 – 25 mm	30 mm	USH10V3F
15 mm		USH15V3F
20 mm		USH20V3F

### STOMOCUR® Safe konvex transparent


Transparent, soft cover one side, soft outlet

Stoma Opening	Plateau Size	Article Number
10 – 43 mm	50 mm	USH10V1FT
10 – 35 mm	40 mm	USH10V2FT
10 – 25 mm	30 mm	USH10V3FT





Available in 3 plateau sizes



### STOMOCUR® Twist & Safe konvex

Window, soft cover both sides, tap valve, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 10 pieces



"Twist and Safe" tap valve

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1H
15 mm		USH15V1H
20 mm		USH20V1H
25 mm		USH25V1H
30 mm		USH30V1H
35 mm		USH35V1H
10 - 35 mm	40 mm	USH10V2H
15 mm		USH15V2H
20 mm		USH20V2H
25 mm		USH25V2H
30 mm		USH30V2H
10 - 25 mm	30 mm	USH10V3H
15 mm		USH15V3H
20 mm		USH20V3H


### STOMOCUR® Twist & Safe konvex transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1HT
10 - 35 mm	40 mm	USH10V2HT
10 - 25 mm	30 mm	USH10V3HT



Available in 3 plateau sizes



### STOMOCUR® Safe medium konvex

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 10 pieces



Soft and comfortable

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1F
15 mm		USH15M1F
20 mm		USH20M1F
25 mm		USH25M1F
30 mm		USH30M1F
35 mm		USH35M1F
10 - 40 mm	45 mm	USH10M2F
15 mm		USH15M2F
20 mm		USH20M2F
25 mm		USH25M2F
10 - 28 mm	35 mm	USH10M3F

### STOMOCUR® Safe medium konvex transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1FT
10 - 40 mm	45 mm	USH10M2FT
10 - 28 mm	35 mm	USH10M3FT

Available  
in 3  
plateau sizes



"Twist and  
Safe" tap  
valve

### STOMOCUR® Twist & Safe medium konvex

Window, soft cover both sides, tap valve,  
adhesive surface 118 mm (Ø), pouch size 130 × 195 mm,  
capacity 420 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
10 – 47 mm	55 mm	USH10M1H
15 mm		USH15M1H
20 mm		USH20M1H
25 mm		USH25M1H
30 mm		USH30M1H
35 mm		USH35M1H
10 – 40 mm	45 mm	USH10M2H
15 mm		USH15M2H
20 mm		USH20M2H
25 mm		USH25M2H
10 – 28 mm	35 mm	USH10M3H

### STOMOCUR® Twist & Safe medium konvex transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Plateau Size	Article Number
10 – 47 mm	55 mm	USH10M1HT
10 – 40 mm	45 mm	USH10M2HT
10 – 28 mm	35 mm	USH10M3HT



Soft and  
comfortable

### STOMOCUR® Safe Clic

Window, soft cover both sides, soft outlet,  
pouch size 130 × 195 mm, capacity 420 ml,  
package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CU45F
57	57 mm	CU57F

### STOMOCUR® Safe Clic transparent

Transparent, soft cover one side, soft outlet

	Locking Ring Size	Article Number
45	45 mm	CU45FT
57	57 mm	CU57FT



"Twist and  
Safe" tap  
valve

### STOMOCUR® Twist & Safe Clic

Window, soft cover both sides, tap valve,  
pouch size 130 × 195 mm, capacity 420 ml,  
package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CU45H
57	57 mm	CU57H

### STOMOCUR® Twist & Safe Clic transparent

Transparent, soft cover one side, tap valve

	Locking Ring Size	Article Number
45	45 mm	CU45HT
57	57 mm	CU57HT

To select the suitable baseplates, please see  
page 29 (STOMOCUR® Clic) and page 30  
(STOMOCUR® Clic konvex).



**STOMOCUR® Safe Microskin**

Window, soft cover both sides, soft outlet, adhesive surface 135 x 100 mm, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
13 - 50 mm	USM13F

**STOMOCUR® Safe Microskin transparent**

Transparent, soft cover one side, soft outlet

Stoma Opening	Article Number
13 - 50 mm	USM13FT



**STOMOCUR® Twist & Safe Microskin**

Window, soft cover both sides, tap valve, adhesive surface 135 x 100 mm, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
13 - 50 mm	USM13H

**STOMOCUR® Twist & Safe Microskin transparent**

Transparent, soft cover one side, tap valve

Stoma Opening	Article Number
13 - 50 mm	USM13HT



**STOMOCUR® Safe gelatinefrei transparent**

Transparent, soft cover one side, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
15 - 60 mm	USG15FT



Accessories



**STOMOCUR® Uro Adapter**

For connecting bed bags to the STOMOCUR® urostomy pouches, package contents 5 pieces

Article Number
USZ1



**STOMOCUR® Drip Cap**

For use with the STOMOCUR® 'Twist and Safe' tap, package contents 5 pieces

Article Number
USZ2



### KONTICUR® Urostomy Bed Bag

Rotary tap drain for one-handed operation, kink-resistant tube, integrated adapter for connection to all FORLIFE urostomy pouches, package contents 5 pieces

Tube Length	Capacity	Article Number
120 cm	2000 ml	KUBB2000



For ostomy patients who do not use products with gelatine or other animal-based ingredients, we offer our pouch assortment with gelatine-free skin protection.

Available for flat one-piece colostomy, ileostomy and urostomy pouches.

#### Features:

- Gentle on the skin
- High tack adhesion
- Thin and flexible
- Hypoallergenic adhesive
- Free of animal components
- Easy to use

#### Benefits:

- Large adhesive area:  
Even when cut to its maximum of 60 mm, there is still enough adhesive to ensure a safe and secure stick
- Large visible release tab offering ease of removal





**STOMOCUR® Select gelatinefrei**

Window, filter, soft cover both sides, adhesive surface 110 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSG20



**STOMOCUR® Protect gelatinefrei**

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 110 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	IPG20

**Soft and discreet closure**



**SPECIAL SUPPLY**

Our comprehensive range of products for specialised care has been developed for medical professionals and includes various types of drainage pouches as well as pouches for paediatric care.

**Paediatric Care:**

- Two options for premature babies up to 900 g and infants up to 5000 g
- Gentle and skin-friendly hydrocolloid
- Soft non-woven cover on one side, transparent back for easy observation
- Closed pouch can be cut open and used as a drainable pouch with a clip or an adhesive tape

**Drainage Protect Pouches:**

- Especially for high output stomata
- Available in flat, convex and medium convex options as one-piece and two-piece pouches
- MPC-Filter (Moisture Protect Composition) for odour control
- High output drain with bung closure
- Hypoallergenic adhesive
- Tapered edge with integral border

**Drainage Pouches:**

- Hypoallergenic adhesive
- Large, flexible skin protection plates combined with flexible, transparent pouches and integrated access port – ideal for promoting healing with easy observation
- High-output drainage ports with soft outlet
- A range of shapes and sizes to meet the patients' needs

STOMOCUR® Safe gelatinefrei transparent - Flat urostomy pouches with gelatine-free skin protection can be found on page 41.







For premature babies up to 900 g

**STOMOCUR® Premature**

Transparent, filter, soft cover one side, adhesive surface 40 mm, pouch size 110 mm x 45 mm, capacity 20 ml, package contents 20 pieces

Stoma Opening	Article Number
3 mm	P3003



**STOMOCUR® Drainage Extra "Access Port"**

Extra Small, opening 10 mm, transparent, package contents 1 piece

Max. Cutability	Access Opening	Article Number
90 x 180 mm	ø 120 mm	SDE2-XSA



**STOMOCUR® Drainage Extra "Access Port"**

Small, opening 10 mm, transparent, package contents 1 piece

Max. Cutability	Access Opening	Article Number
240 x 155 mm	240 x 155 mm	SDE2-SA



**STOMOCUR® Drainage Extra "Access Port"**

Medium, opening 10 mm, transparent, package contents 1 piece

Max. Cutability	Access Opening	Article Number
85 x 250 mm	85 x 250 mm	SDE2-MA



**STOMOCUR® Drainage Extra "Access Port"**

Large, opening 10 mm, transparent, package contents 1 piece

Max. Cutability	Access Opening	Article Number
85 x 250 mm	160 x 200 mm	SDE2-LA



For infants up to 5000 g

**STOMOCUR® Premature Combi**

Transparent, soft cover one side, adhesive surface 40 mm, pouch size 210 mm x 45 mm, capacity 90 ml, package contents 30 pieces

Stoma Opening	Article Number
3 mm	P3203



With 2 drainage outlets



**STOMOCUR® Drainage Extra "Access Port"**

Extra Large, opening 10 mm, foil: transparent, package contents 1 piece

Max. Cutability	Access Opening	Article Number
230 x 380 mm	ø 120 mm	SDE2-XXLA

Soft and flexible



**STOMOCUR® Drainage Protect**

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Article Number
20 - 60 mm	DP5020

**STOMOCUR® Drainage small - large**

Transparent with printed scale, backstop, exhaust valve with adapter, package contents 30 pieces

Stoma Opening	Capacity	Article Number
3 mm	50 ml	DP2503XS
5 mm	150 ml	DP4005S
5 mm	300 ml	DP4005



**STOMOCUR® Drainage Protect konvex**

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, convex height 7 mm, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	DP20V1
15 - 35 mm	40 mm	DP15V2
10 - 25 mm	30 mm	DP10V3

Available in 3 plateau sizes



**STOMOCUR® Drainage Protect medium konvex**

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, convex height 3,5 mm, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	DP20M1
15 - 40 mm	45 mm	DP15M2
10 - 28 mm	35 mm	DP10M3

Available in 3 plateau sizes





### STOMOCUR® Drainage Protect Bed Bag

Large-lumen and kink-resistant tube, integrated adapter for connection to all FORLIFE Drainage Protect pouches, package contents 5 pieces

Tube Length	Capacity	Article Number
130 mm	2.000 ml	DPB2000

#### Compatible with:

STOMOCUR® Drainage Protect  
 STOMOCUR® Drainage Protect konvex  
 STOMOCUR® Drainage Protect medium konvex  
 STOMOCUR® Drainage Protect Clic

Stoma Opening	Plateau Size	Article Number
20 – 60 mm	-	DP5020
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3
20 – 47 mm	55 mm	DP20M1
15 – 40 mm	45 mm	DP15M2
10 – 28 mm	35 mm	DP10M3
Locking ring: 70 mm	-	CDP70T




The comprehensive FORLIFE range of ostomy care products and accessories provides additional safety and reassurance.

- Hydrocolloid skin protection plates
- Skin protection rings
- Fixing strips
- Skin protection paste
- Skin protection powder
- Nourishing lotion
- Adhesive remover spray
- Ostomy belts
- Superabsorbent gelling agents
- Urostomy adapters
- Urostomy drip cap
- Scissors



## Accessories



### STOMOCUR® Skin Protection Plates

With self-adhesive hydrocolloid skin protection compound, material thickness 0,6 mm, package contents 5 pieces

Size	Article Number
100 x 100 mm	HP101006
150 x 150 mm	HP151506
200 x 200 mm	HP202006
370 x 220 mm	HP372206

### STOMOCUR® Skin Protection Plates mouldable

With self-adhesive hydrocolloid skin protection compound, can be individually cut and shaped, material thickness 2,4 mm, package contents 1 piece

Size	Article Number
220 x 220 mm	HP222224

### STOMOCUR® Skin Protection Rings

With hydrocolloid skin protection material to seal and to protect the skin around the stoma, available in two different outer diameters, material thickness 2,4 mm, package contents 30 pieces

Outer Ø	Size	Article Number
48	10 mm	48101
48	15 mm	48151
48	20 mm	48201
63	25 mm	63251
63	30 mm	63301
63	35 mm	63351
63	40 mm	63401



## Accessories



### STOMOCUR® Fixing Strips

Skin-friendly and highly elastic, extends the wearing time, prevents the baseplate from detaching and rolling up, more safety and comfort, package contents 20 pieces

Article Number
FS20



### STOMOCUR® Skin Protection Paste

For external use on the skin around the stoma, fills unevenness between the stoma and the skin protection plate, available with and without alcohol, package contents 1 piece

Article Number	Contents
CP2050 (with alcohol)	56 g
PASTEOA (without alcohol)	56 g



### STOMOCUR® Skin Protection Powder

The powder absorbs the moisture of the skin around the stoma, ensuring better adhesion of the stoma supply, package contents 1 piece

Article Number	Contents
HP30	30 g





## Accessories



### STOMOCUR® Safety

The skin care lotion is applied very thinly to the skin around the stoma to prevent skin irritation, package contents 1 piece

Article Number	Contents
SSAFETY100	100 ml



### STOMOCUR® EmplaSectal

For easy and skin-friendly removal of the skin protection plates, handy as a spray, package contents 1 piece

Article Number	Contents
EMPLASECTAL	50 ml



### STOMOCUR® Ostomy Belt

Ostomy belt for additional safety, available in two lengths, adjustable, width 3 cm, package contents 1 piece

Article Number	Length
AB01	100 cm
AB02	120 cm



### STOMOCUR® SGX

Fast-acting, super-absorbent gelling agent, easy to use, helpful for ileostomists and colostomists with recipients with more fluid output, reduces pouch swelling, pouch noise and odour, ideal for travel, package contents 30 pieces

Article Number
SGX30



## Accessories



### STOMOCUR® Uro Adapter

For connecting bed bags to the STOMOCUR® urostomy pouches, package contents 5 pieces

Article Number
USZ1



### STOMOCUR® Uro Drip Cap

To attach the STOMOCUR® 'Twist and Safe' valve, package contents 5 pieces

Article Number
USZ2



### STOMOCUR® Scissors

Scissors with rounded blades for cutting the baseplate openings, package contents 1 piece

Article Number
SSchere



Table of contents by item number

Article Number	Detail	Page
48101	10 mm	52
48151	15 mm	52
48201	20 mm	52
63251	25 mm	52
63301	30 mm	52
63351	35 mm	52
63401	40 mm	52
AB01		54
AB02		54
CA4515V2	15 mm - 35 mm	30
CA4515V3	15 mm - 25 mm	30
CA5715V1	15 mm - 43 mm	30
CA5715V2	15 mm - 35 mm	30
CA5715V3	15 mm - 25 mm	30
CA4515	15 mm - 35 mm	29
CA5715	15 mm - 40 mm	29
CAA20V1B	20 mm - 43 mm	10
CAA15V2B	15 mm - 35 mm	10
CAA10V3B	10 mm - 25 mm	10
CAA10M3B	10 mm - 28 mm	12
CAA15M2B	15 mm - 40 mm	12
CAA20M1B	20 mm - 47 mm	12
CAH10V3B	10 mm - 25 mm	9
CAH15V2B	15 mm - 35 mm	9
CAH20V1B	20 mm - 43 mm	9
CAH10M3B	10 mm - 28 mm	11
CAH15M2B	15 mm - 40 mm	11
CAH20M1B	20 mm - 47 mm	11
CC45	45 mm	31
CC57	57 mm	31
CC70	70 mm	31
CDP70T	70 mm	32 & 50
CI45	45 mm	31
CI45S	45 mm	32
CI57	57 mm	31
CI70	70 mm	31
CK4515	15 mm - 35 mm	29
CK4515V2	15 mm - 35 mm	30
CK4515V3	15 mm - 25 mm	30
CK4520	20 mm	29
CK4520V2	20 mm	30
CK4525	25 mm	29
CK4525V2	25 mm	30
CK4530	30 mm	29
CK4530V2	30 mm	30
CK5715	15 mm - 40 mm	29
CK5715V1	15 mm - 43 mm	30
CK5715V2	15 mm - 35 mm	30
CK5715V3	15 mm - 25 mm	30
CK5720V1	20 mm	30
CK5725V1	25 mm	30

Table of contents by item number

Article Number	Detail	Page
CK5730V1	30 mm	30
CK5735	35 mm	29
CK5735V1	35 mm	30
CK5740	40 mm	29
CK7015	15 mm - 50 mm	29
CK7045	45 mm	29
CK7050	50 mm	29
CP2050		53
CSA10V3	10 mm - 25 mm	20
CSA15V2	15 mm - 35 mm	20
CSA15V3	15 mm	20
CSA20	20 mm - 60 mm	26
CSA20S	20 mm - 60 mm	26
CSA20V1	20 mm - 43 mm	20
CSA20V2	20 mm	20
CSA20V3	20 mm	20
CSA25	25 mm	26
CSA25S	25 mm	26
CSA25V1	25 mm	20
CSA25V2	25 mm	20
CSA30	30 mm	26
CSA30S	30 mm	26
CSA30V1	30 mm	20
CSA30V2	30 mm	20
CSA35	35 mm	26
CSA35S	35 mm	26
CSA35V1	35 mm	20
CSA40	40 mm	26
CSA40S	40 mm	26
CSA40V1	40 mm	20
CSC20	20 mm - 50 mm	17
CSC20SNF	20 mm - 40 mm	17
CSC25NF	25 mm - 50 mm	17
CSC30LNF	30 mm - 60 mm	17
CSG20	20 mm - 60 mm	44
CSH10M3	10 mm - 28 mm	23
CSH10V3	10 mm - 25 mm	19
CSH10XS	10 mm - 50 mm	15
CSH10XST	10 mm - 50 mm	15
CSH15M2	15 mm - 40 mm	23
CSH15V2	15 mm - 35 mm	19
CSH15V3	15 mm	19
CSH20	20 mm - 60 mm	14
CSH20L	20 mm - 90 mm	14
CSH20M1	20 mm - 47 mm	23
CSH20M2	20 mm	23
CSH20S	20 mm - 60 mm	15
CSH20V0	20 mm - 57 mm	19
CSH20V1	20 mm - 43 mm	19
CSH20V2	20 mm	19
CSH20V3	20 mm	19



Table of contents by item number

Article Number	Detail	Page
CSH25	25 mm	14
CSH25L	25 mm	14
CSH25M1	25 mm	23
CSH25M2	25 mm	23
CSH25S	25 mm	15
CSH25V1	25 mm	19
CSH25V2	25 mm	19
CSH30	30 mm	14
CSH30L	30 mm	14
CSH30M1	30 mm	23
CSH30S	30 mm	15
CSH30V1	30 mm	19
CSH30V2	30 mm	19
CSH35	35 mm	14
CSH35L	35 mm	14
CSH35M1	35 mm	23
CSH35S	35 mm	15
CSH35V1	35 mm	19
CSH40	40 mm	14
CSH40L	40 mm	14
CSH40M1	40 mm	23
CSH40S	40 mm	15
CSH40V1	40 mm	19
CU45F	45 mm	39
CU45FT	45 mm	39
CU45H	45 mm	39
CU45HT	45 mm	39
CU57F	57 mm	39
CU57FT	57 mm	39
CU57H	57 mm	39
CU57HT	57 mm	39
DP10M3	10 mm - 28 mm	24, 49 & 50
DP10V3	10 mm - 25 mm	21, 49 & 50
DP15M2	15 mm - 40 mm	24, 49 & 50
DP15V2	15 mm - 35 mm	21, 49 & 50
DP20M1	20 mm - 47 mm	24, 49 & 50
DP20V1	20 mm - 43 mm	21, 49 & 50
DP2503XS	3 mm	48
DP4005	5 mm	48
DP4005S	5 mm	48
DP5020	20 mm - 60 mm	17, 49 & 50
DPB2000		50
EMPLASECTAL		54
FS20		53
HP101006	100 mm x 100 mm	52
HP151506	150 mm x 150 mm	52
HP202006	200 mm x 200 mm	52
HP222224	220 mm x 220 mm	52
HP30		53
HP372206	370 mm x 220 mm	52
IAA10M3B	10 mm - 28 mm	12

Table of contents by item number

Article Number	Detail	Page
IAA15M2B	15 mm - 40 mm	12
IAA15M1B	15 mm - 47 mm	12
IAA10V3B	10 mm - 25 mm	10
IAA15V1B	15 mm - 43 mm	10
IAA15V2B	15 mm - 35 mm	10
IAH10M3B	10 mm - 28 mm	11
IAH15M1B	15 mm - 47 mm	11
IAH15M2B	15 mm - 40 mm	11
IPA10V3	10 mm - 25 mm	20
IPA15	15 mm - 60 mm	27
IPA15V1	15 mm - 43 mm	20
IPA15V2	15 mm - 35 mm	20
IPA20	20 mm	27
IPA20V1	20 mm	20
IPA25	25 mm	27
IPA25V1	25 mm	20
IPA30	30 mm	27
IPA30V1	30 mm	20
IPA35	35 mm	27
IPA35V1	35 mm	20
IPA40	40 mm	27
IPA40V1	40 mm	20
IPG20	20 mm - 60 mm	44
IPH10M3	10 mm - 28 mm	23
IPH10M3T	10 mm - 28 mm	23
IPH10V3	10 mm - 25 mm	19
IPH15	15 mm - 60 mm	16
IPH15M1	15 mm - 47 mm	23
IPH15M1T	15 mm - 47 mm	23
IPH15M2	15 mm - 40 mm	23
IPH15M2T	15 mm - 40 mm	23
IPH15T	15 mm - 60 mm	16
IPH15S	15 mm - 60 mm	16
IPH15V1	15 mm - 43 mm	19
IPH15V2	15 mm - 35 mm	19
IPH15V3	15 mm	19
IPH20	20 mm	16
IPH20M1	20 mm	23
IPH20M2	20 mm	23
IPH20V0	20 mm - 57 mm	19
IPH20V1	20 mm	19
IPH20V2	20 mm	19
IPH20V3	20 mm	19
IPH2070	20 mm - 70 mm	16
IPH25	25 mm	16
IPH25M1	25 mm	23
IPH25M2	25 mm	23
IPH25V1	25 mm	19
IPH25V2	25 mm	19
IPH30	30 mm	16
IPH30M1	30 mm	23

Table of contents by item number

Article Number	Detail	Page
IPH30V1	30 mm	19
IPH30V2	30 mm	19
IPH35	35 mm	16
IPH35M1	35 mm	23
IPH35V1	35 mm	19
IPH40	40 mm	16
IPH40M1	40 mm	23
IPH40V1	40 mm	19
KUBB2000		42
P3003	3 mm	46
P3203	3 mm	46
PASTE0A		53
SDE2-LA	160 mm × 200 mm	47
SDE2-MA	85 mm × 250 mm	47
SDE2-SA	240 mm × 155 mm	47
SDE2-XSA	90 mm × 180 mm	47
SDE2-XXLA	230 mm × 380 mm	48
SGX30		54
SSAFETY100		54
SSchere		55
USG15FT	15 mm – 60 mm	41
USH10F	10 mm – 60 mm	34
USH10FT	10 mm – 60 mm	34
USH10H	15 mm – 60 mm	34
USH10HT	10 mm – 60 mm	34
USH10M1F	10 mm – 47 mm	37
USH10M1FT	10 mm – 47 mm	37
USH10M1H	10 mm – 47 mm	38
USH10M1HT	10 mm – 47 mm	38
USH10M2F	10 mm – 40 mm	37
USH10M2FT	10 mm – 40 mm	37
USH10M2H	10 mm – 40 mm	38
USH10M2HT	10 mm – 40 mm	38
USH10M3F	10 mm – 28 mm	37
USH10M3FT	10 mm – 28 mm	37
USH10M3H	10 mm – 28 mm	38
USH10M3HT	10 mm – 28 mm	38
USH10V1H	10 mm – 43 mm	36
USH10V1HT	10 mm – 43 mm	36
USH10V1F	10 mm – 43 mm	35
USH10V1FT	10 mm – 43 mm	35
USH10V2F	10 mm – 35 mm	35
USH10V2FT	10 mm – 35 mm	35
USH10V2H	10 mm – 35 mm	36
USH10V2HT	10 mm – 35 mm	36
USH10V3F	10 mm – 25 mm	35
USH10V3FT	10 mm – 25 mm	35
USH10V3H	10 mm – 25 mm	36
USH10V3HT	10 mm – 25 mm	36
USH15F	15 mm	34
USH15H	15 mm	34

Table of contents by item number

Article Number	Detail	Page
USH15M1F	15 mm	37
USH15M1H	15 mm	38
USH15M2F	15 mm	37
USH15M2H	15 mm	38
USH15V1F	15 mm	35
USH15V1H	15 mm	36
USH15V2F	15 mm	35
USH15V2H	15 mm	36
USH15V3F	15 mm	35
USH15V3H	15 mm	36
USH20F	20 mm	34
USH20H	20 mm	34
USH20M1F	20 mm	37
USH20M1H	20 mm	38
USH20M2F	20 mm	37
USH20M2H	20 mm	38
USH20V1F	20 mm	35
USH20V1H	20 mm	36
USH20V2F	20 mm	35
USH20V3F	20 mm	35
USH20V2H	20 mm	36
USH20V3H	20 mm	36
USH25F	25 mm	34
USH25H	25 mm	34
USH25M1F	25 mm	37
USH25M1H	25 mm	38
USH25M2F	25 mm	37
USH25M2H	25 mm	38
USH25V1F	25 mm	35
USH25V1H	25 mm	36
USH25V2F	25 mm	35
USH25V2H	25 mm	36
USH30F	30 mm	34
USH30H	30 mm	34
USH30M1F	30 mm	37
USH30M1H	30 mm	38
USH30V1F	30 mm	35
USH30V1H	30 mm	36
USH30V2F	30 mm	35
USH30V2H	30 mm	36
USH35M1F	35 mm	37
USH35M1H	35 mm	38
USH35V1F	35 mm	35
USH35V1H	35 mm	36
USM13F	13 mm – 50 mm	40
USM13FT	13 mm – 50 mm	40
USM13H	13 mm – 50 mm	40
USM13HT	13 mm – 50 mm	40
USZ1		41 & 55
USZ2		41 & 55



For Life  
Produktions- und Vertriebsgesellschaft  
für Heil- und Hilfsmittel mbH

Am Studio 16  
12489 Berlin, Germany

[mail@forlife.info](mailto:mail@forlife.info)  
[www.forlife.info](http://www.forlife.info)



**FORLIFE**  
BERLIN

*Empowering yourself.*