



We are expanding our two-piece **STOMOCUR® Clic System with two new special features!**

- Now available in our skin protecting Alginate range.
- The unique, improved locking ring is even easier to handle from behind to reduce pressure when attaching to the pouch.



Hypoallergenic

Healing

Absorbent

Antibacterial

Haemostatic



Improved handling of the locking ring system

- A simple and stable connection between the baseplate and the pouch.
- An audible click confirms correct attachment, resulting in increased safety.

Why choose our two-piece system?

- Flexible and easy to connect system
- Available in flat and convex baseplates
- Soft non-woven pouch cover for additional comfort
- Wide ranges of different pouch types
- Three different locking ring sizes (45 mm, 57 mm or 70 mm)

Two-Piece Supply | Clic System | Baseplates



STOMOCUR® Clic & Clic Alginate

Adhesive surface 118 mm (Ø), package contents 5 pieces

Locking ring	Stoma opening	Article number
45	15 - 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
NEW	15 - 35 mm	CA4515
57	15 - 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
NEW	15 - 40 mm	CA5715
70	15 - 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050

Alginat Hydrocolloid



STOMOCUR® Clic konvex & Clic Alginate konvex


Adhesive surface 118 mm (Ø), package contents 3 pieces

Locking ring	Stoma opening	Plateau size	Article number
45	15 - 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
NEW	15 - 35 mm		CA4515V2
NEW	15 - 25 mm	30 mm	CK4515V3
	15 - 25 mm	30 mm	CA4515V3
57	15 - 43 mm	50 mm	CK5715V1
	20 mm		CK5720V1
	25 mm		CK5725V1
	30 mm		CK5730V1
	35 mm		CK5735V1
NEW	15 - 43 mm		CA5715V1
NEW	15 - 35 mm	40 mm	CK5715V2
	15 - 35 mm	40 mm	CA5715V2
NEW	15 - 25 mm	30 mm	CK5715V3
	15 - 25 mm	30 mm	CA5715V3

Two-Piece Supply | Clic System | Pouches

Package contents 30 pieces

	Article name	Pouch type	Article number
45	STOMOCUR® Select Clic	Colostomy	CC45
	STOMOCUR® Protect Clic	Ileostomy	CI45
	STOMOCUR® Protect Clic small	Ileostomy	CI45S
	STOMOCUR® Safe Clic	Urostomy	CU45F
	STOMOCUR® Safe Clic transparent	Urostomy	CU45FT
	STOMOCUR® Twist & Safe Clic	Urostomy	CU45H
	STOMOCUR® Twist & Safe Clic transparent	Urostomy	CU45HT
57	STOMOCUR® Select Clic	Colostomy	CC57
	STOMOCUR® Protect Clic	Ileostomy	CI57
	STOMOCUR® Safe Clic	Urostomy	CU57F
	STOMOCUR® Safe Clic transparent	Urostomy	CU57FT
	STOMOCUR® Twist & Safe Clic	Urostomy	CU57H
	STOMOCUR® Twist & Safe Clic transparent	Urostomy	CU57HT
70	STOMOCUR® Select Clic	Colostomy	CC70
	STOMOCUR® Protect Clic	Ileostomy	CI70
	STOMOCUR® Drainage Protect Clic transparent	Drainage	CDP70T



Date: June 2023