

We are expanding our one-piece ileostomy range

STOMOCUR® Protect

- Our proven one-piece flat ileostomy pouch, STOMOCUR® Protect, is now available in cut-to-fit up to 70 mm!
- New STOMOCUR® Protect IPH2070, ideal for very large stomas or for post-op supplies.

Our advantages:

Hypoallergenic continuous hydrocolloid skin protection

Good skin tolerance

High degree of flexibility ensures a comfortable hold

High absorption capacity

Patented MPC filter technology with additional liquid protection layer

Unique, easy-to-use 'Safe and Easy' outlet

Soft non-woven cover on both sides

Better visibility due to 'Window' technology

cut-to-fit
70mm



STOMOCUR® Protect

Adhesive surface 110 mm Ø, capacity 530 ml, pack contains 30 pieces

Stoma opening	Article number
15 – 60 mm	IPH15
20 mm	IPH20
25 mm	IPH25
30 mm	IPH30
35 mm	IPH35
40 mm	IPH40

NEW! 20 – 70 mm

IPH2070

