



STOMOCUR[®]

Product Guide 2024
Stomacare

Introduction

FORLIFE has stood for high-quality products for over 30 years, enabling economic activities for medical retailers, clinics and pharmacies while maintaining the highest quality standards – Made in Germany.

We have deliberately chosen to conduct research, development and production exclusively in Germany. As well as trusting in the extensive experience of our employees, the location in Germany is also part of our corporate philosophy. We are delighted that we have successfully brought together long-standing FORLIFE employees with younger talents, combining proven technologies with innovation.

Progress is an elementary component of our corporate success. In addition to the expertise of our employees, we rely on state-of-the-art production machinery. The use of cutting-edge technologies and efficient production methods ensure the long-term economic production of high-quality and innovative products in Germany.

Our goal is to provide ideal care for the stoma, allowing ostomates to get on with their lives in comfort and security. This is the objective of our products and our daily work.

Your FORLIFE Team



Always well taken care of with STOMOCUR®

Our STOMOCUR® pouches are manufactured in Germany and characterized by our innovative and unique features. We use skin protection materials made from natural ingredients that promote skin health and reduce irritation and inflammation.

Our pouches are comfortable and secure to wear and provide confidence and comfort for the user. Ultimately ensuring healthy skin around the stoma.

All our products are made from a very soft, quick-drying non-woven fabric. Our pouches contain Polyethylene (PE), which does not stick to clothing and is gentle on the skin to avoid skin reactions. Our skin protection materials are composed of natural ingredients that both promote skin health and reduce irritation and inflammation.

Many of our pouches are equipped with a 'Window', allowing patients and medical professionals to observe the stoma without removing the pouch.

To meet individual needs, STOMOCUR® pouches are available as a one-piece system, where the baseplate is connected to the pouch and as a two-piece 'Clic' system, where the pouch can be removed and replaced without replacing the baseplate.

Our Colostomy 'Select' range has a pouch shape resembling a symmetrical figure-of-eight. This shape fits on either side of the body and is designed for maximum capacity with minimum bulk, allowing the product to be worn for longer without being noticeable under clothing.

The unique activated carbon filter is an integral part of all our colostomy pouches, ensuring odour-free use for the user.

Our Ileostomy drainable pouches, 'Protect', are equipped with our patented MPC filter technology (Moisture Protect Composition) - an additional protective layer that prevents liquid absorption, keeps the user odour-free and allows unwanted gases to escape.

All drainable pouches feature our unique, easy-to-use 'Safe and Easy' outlet.


Our Urostomy pouches are available with soft bung closure, the 'Safe' range, as well as our tap closure, 'Twist & Safe' both providing choice, comfort and confidence. Night drain adapters are included with all urostomy pouches and 'Twist & Safe' pouches also come with a drip cap to prevent leaks.

Further details on the respective product ranges can be found on the following pages of our STOMOCUR® product guide.



Convex care

For ultimate safety and comfort, we offer a range of convex care supplies. STOMOCUR® offers a total of four different convexities in two different heights. This range of plateau sizes gives the user maximum choice.



Contents



HYDROKOLLOID

Hydrocolloid

Soft, flexible, hypoallergenic skin protection.

Page 8



COMBIFLANGE

Combiflange Hydrocolloid

Contoured hydrocolloid for excellent adhesion and assured safety against leaks.

Page 8



Convex

7 mm convexity. Available with contoured Combiflange hydrocolloid or Alginate skin protection. Also suitable for drainage pouches with Combiflange hydrocolloid skin protection.

Page 13



Urostomy

Available with tap and soft bung closure options with different skin protection materials.

Page 28



Microskin

The thinnest and most flexible skin protection.

Page 35



gelatinefrei

Gelatine Free

Using a hydrocolloid skin protection without gelatine or other animal components. Gentle on the skin.

Page 37



Medium Convex

3,5 mm convexity. Available with contoured Combiflange hydrocolloid. Also suitable for drainage pouches with Combiflange hydrocolloid skin protection.

Page 17



Alginate

An alginate-containing skin protection. Promotes skin integrity, absorbs moisture, and aids in healing damaged skin.

Page 20



Clic System

A versatile, flexible and discreet two-piece system.

Page 23



Special Supply

Wound and fistula drainage supply in various sizes.

Page 39



Accessories

A comprehensive range of accessories.

Page 44



HYDROKOLLOID



COMBIFLANGE

Hydrocolloid and Combiflange

FORLIFE's Hydrocolloid is a flexible, hypoallergenic, and skin-friendly skin protection, ideal for stoma care.

How does Hydrocolloid work?

Hydrocolloid adheres to the skin through the absorption of heat and moisture. This moisture absorption causes the baseplate to swell slightly, ensuring good adhesion to the skin and preventing leaks and skin irritation.

Features:

- Hypoallergenic skin protection
- Good skin compatibility
- Comfortable hold due to maximum flexibility
- Easy application

Combiflange:

Baseplates with Combiflange hydrocolloid are made up of two interconnected hydrocolloid layers that are thicker on the inside, directly on the stoma, and become thinner towards the outside. This double layer ensures a quick and secure hold that lasts longer and has a high absorption capacity. The Combiflange skin protection is available in a variety of our convex STOMOCUR® products, flat urostomy and drainage pouches as well as two-piece convex and flat baseplates.



One-Piece Supply



STOMOCUR® Select large

Window, filter, soft cover both sides, adhesive surface 135 × 135 mm, pouch size 160 × 220 mm, capacity 470 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 90 mm	CSH20L
25 mm	CSH25L
30 mm	CSH30L
35 mm	CSH35L
40 mm	CSH40L



STOMOCUR® Select

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 145 × 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSH20
25 mm	CSH25
30 mm	CSH30
35 mm	CSH35
40 mm	CSH40



One-Piece Supply

**STOMOCUR® Select small**

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 130 × 190 mm, capacity 250 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSH20S
25 mm	CSH25S
30 mm	CSH30S
35 mm	CSH35S
40 mm	CSH40S

**STOMOCUR® Select extra small**

Window, filter, soft cover both sides, adhesive surface 80 mm (Ø), pouch size 100 × 160 mm, capacity 140 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 50 mm	CSH10XS

STOMOCUR® Select extra small transparent

Transparent, filter, soft cover one side

Stoma Opening	Article Number
10 – 50 mm	CSH10XST

Also suitable for child care

**STOMOCUR® Select Cap**

Adhesive surface 60 mm (Ø), package contents 30 pieces

Stoma Opening	Article Number
20 – 50 mm	CSC20

After Irrigation



One-Piece Supply

Easy to apply

**STOMOCUR® Protect**

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 110 mm (Ø), pouch size 150 × 220 mm, capacity 530 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 70 mm	IPH2070
15 – 60 mm	IPH15
20 mm	IPH20
25 mm	IPH25
30 mm	IPH30
35 mm	IPH35
40 mm	IPH40

STOMOCUR® Protect transparent

Transparent, filter, soft cover one side, 'Safe and Easy' closure

Stoma Opening	Article Number
15 – 60 mm	IPH15T

**STOMOCUR® Protect small**

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 110 mm (Ø), pouch size 150 × 190 mm, capacity 400 ml, package contents 30 pieces

Stoma Opening	Article Number
15 – 60 mm	IPH15S

Also suitable for child care



Special Supply

**STOMOCUR® Drainage Protect**

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, package contents 10 pieces


Stoma Opening	Article Number
20 – 60 mm	DP5020





Our convex ostomy pouches are the flexible choice for maximum security and additional confidence. They adapt to the contours of your body, ensuring more comfort during wear. Combiflange hydrocolloid and alginate options are available in our Konvex range.

Features:

- 7 mm convexity
- Maintains its convex shape during wear
- Easy-to-use 'Safe and Easy' outlet for drainable pouches
- Comprehensive range of cut-to-fit or pre-cut sizes
- High adhesion for advanced performance
- Skin-friendly, hypoallergenic skin protection
- Efficient odour elimination through the integrated filter system
- Four convexities
- Integrated belt attachment

Convex supply:



STOMOCUR® Select konvex

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces



Stoma Opening	Plateau Size	Article Number
20 – 57 mm	60 mm	CSH20V0
20 – 43 mm	50 mm	CSH20V1
25 mm		CSH25V1
30 mm		CSH30V1
35 mm		CSH35V1
40 mm		CSH40V1
15 – 35 mm	40 mm	CSH15V2
20 mm		CSH20V2
25 mm		CSH25V2
30 mm		CSH30V2
10 – 25 mm	30 mm	CSH10V3
15 mm		CSH15V3
20 mm		CSH20V3

STOMOCUR® Select Alginat konvex

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces



For sensitive skin

Stoma Opening	Plateau Size	Article Number
20 – 43 mm	50 mm	CSA20V1
25 mm		CSA25V1
30 mm		CSA30V1
35 mm		CSA35V1
40 mm		CSA40V1
15 – 35 mm	40 mm	CSA15V2
20 mm		CSA20V2
25 mm		CSA25V2
30 mm		CSA30V2
10 – 25 mm	30 mm	CSA10V3
15 mm		CSA15V3
20 mm		CSA20V3

STOMOCUR® Protect konvex

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 10 pieces



Easy to apply

Stoma Opening	Plateau Size	Article Number
20 – 57 mm	60 mm	IPH20V0
15 – 43 mm	50 mm	IPH15V1
20 mm		IPH20V1
25 mm		IPH25V1
30 mm		IPH30V1
35 mm		IPH35V1
40 mm	IPH40V1	
15 – 35 mm	40 mm	IPH15V2
20 mm		IPH20V2
25 mm		IPH25V2
30 mm		IPH30V2
10 – 25 mm	30 mm	IPH10V3
15 mm		IPH15V3
20 mm		IPH20V3

STOMOCUR® Protect Alginat konvex

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 30 ml, package contents 10 pieces




For sensitive skin

Stoma Opening	Plateau Size	Article Number
15 – 43 mm	50 mm	IPA15V1
20 mm		IPA20V1
25 mm		IPA25V1
30 mm		IPA30V1
35 mm		IPA35V1
40 mm	IPA40V1	
15 – 35 mm	40 mm	IPA15V2
10 – 25 mm	30 mm	IPA10V3



Special Supply

Available in 3 plateau sizes



Soft and flexible

STOMOCUR® Drainage Protect konvex

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3




Our Medium Konvex range offers the same outstanding flexibility and comfort as Konvex but with a lower convexity of 3,5 mm. Combiflange hydrocolloid is available as convex skin protection.

Features:

- Low profile: 3,5 mm convexity
- Comprehensive range of cut-to-fit or pre-cut sizes
- High adhesion for advanced performance
- Skin-friendly, hypoallergenic skin protection
- Efficient odor elimination through the integrated filter system
- Easy-to-use 'Safe and Easy' outlet for drainable pouches
- Soft and flexible
- Three convexities
- Integrated belt attachment

Medium convex supply:





STOMOCUR® Select medium konvex

Window, filter, soft cover both sides, adhesive surface 118 mm (Ø), pouch size 145 × 205 mm, capacity 380 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	CSH20M1
25 mm		CSH25M1
30 mm		CSH30M1
35 mm		CSH35M1
40 mm		CSH40M1
15 - 40 mm	45 mm	CSH15M2
20 mm		CSH20M2
25 mm		CSH25M2
10 - 28 mm	35 mm	CSH10M3



Soft and flexible

STOMOCUR® Drainage Protect medium konvex

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 × 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	DP20M1
15 - 40 mm	45 mm	DP15M2
10 - 28 mm	35 mm	DP10M3



STOMOCUR® Protect medium konvex

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 118 mm (Ø), pouch size 150 × 220 mm, capacity 530 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IPH15M1
20 mm		IPH20M1
25 mm		IPH25M1
30 mm		IPH30M1
35 mm		IPH35M1
40 mm	IPH40M1	
15 - 40 mm	45 mm	IPH15M2
20 mm		IPH20M2
25 mm		IPH25M2
10 - 28 mm	35 mm	IPH10M3

STOMOCUR® Protect medium konvex transparent

Transparent, filter, soft cover one side, 'Safe and Easy' closure

Stoma Opening	Plateau Size	Article Number
15 - 47 mm	55 mm	IPH15M1T
15 - 40 mm	45 mm	IPH15M2T
10 - 28 mm	35 mm	IPH10M3T

Easy to apply

Bioactive skin healing

Our alginate baseplate is one of our top skin protection products, offering both preventive and healing properties. Available in a variety of sizes for both drainable and closed pouches. Alginate has been used in woundcare for many years.

Features:

- Healing
- Absorbent
- Antibacterial
- Haemostatic
- Used without any accessories to cost-effectively treat problematic skin

How does Alginate work?

We exclusively use natural alginate found in brown seaweed, which can absorb up to 20 times its weight in liquid. The alginate in our skin protection absorbs moisture from the skin, making it ideal for healing damaged skin and keeping it healthy.

Applications for use are:

- Skin prone to allergies
- Sore, excoriated peri-stomal skin
- Moist bleeding, peri-stomal skin
- Ulcerated peri-stomal skin
- Muco-cutaneous separation
- Stomas which produce excessive mucous
- Infected peri-stomal suture sites
- Alginate adhesive can readily be used on non-problematic skin to prevent complications occurring

FORLIFE Alginate skin protection contains not just a trace element but 15% Alginate, this is significant enough to promote healing.

One-Piece Supply



For sensitive skin

STOMOCUR® Select Alginat

Window, filter, soft cover both sides, adhesive surface 110 × 110mm, pouch size 145 × 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSA20
25 mm	CSA25
30 mm	CSA30
35 mm	CSA35
40 mm	CSA40



For sensitive skin

STOMOCUR® Select Alginat small

Window, filter, soft cover both sides, adhesive surface 110 × 110 mm, pouch size 130 × 190 mm, capacity 250 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSA20S
25 mm	CSA25S
30 mm	CSA30S
35 mm	CSA35S
40 mm	CSA40S

STOMOCUR® Select Alginat konvex - Convex colostomy pouches with Alginate skin protection can be found on page 15.



One-Piece Supply



STOMOCUR® Protect Alginat

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 110 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces



Stoma Opening	Article Number
15 - 60 mm	IPA15
20 mm	IPA20
25 mm	IPA25
30 mm	IPA30
35 mm	IPA35
40 mm	IPA40



For those who prefer a two-piece ostomy pouch system, we offer our popular Clic System for closed, drainable and drainage pouches, as well as for urostomy care.

With three different ring sizes (45 mm, 57 mm or 70 mm) and with flat or convex baseplates, the Clic System is flexible and adapts to various needs. Combiflange hydrocolloid and alginate are available as skin protection.

Features:

- Versatile, flexible and easy-to-connect system
- Unique spacing between the ring and the baseplate allows the user to attach the baseplate without applying pressure
- Provides a simple, stable connection between the baseplate and pouch
- Audible 'click' when attaching the pouch confirms that it has been attached securely for safety against leaks
- Soft non-woven cover for additional comfort
- Wide range of pouches for both types of baseplates

STOMOCUR® Protect Alginat konvex - Convex ileostomy pouches with Alginate skin protection can be found on page 15.

Two-Piece Supply | Clic System | Baseplates



STOMOCUR® Clic

Plate size 118 mm (Ø), package contents 5 pieces

	Stoma Opening	Article Number
45	15 – 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
57	15 – 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
70	15 – 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050

Two-Piece Supply | Clic System | Baseplates



STOMOCUR® Clic konvex

Plate size 118 mm (Ø), package contents 3 pieces

	Stoma Opening	Plateau Size	Article Number
45	15 – 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
57	15 – 25 mm	30 mm	CK4515V3
	15 – 43 mm	50 mm	CK5715V1
	20 mm		CK5720V1
25 mm	CK5725V1		
30 mm	CK5730V1		
70	35 mm	CK5735V1	
	15 – 35 mm	40 mm	CK5715V2
	15 – 25 mm	30 mm	CK5715V3

Available
in 2
plateau sizes



STOMOCUR® Clic Alginat

Plate size 118 mm (Ø), package contents 5 pieces

	Stoma Opening	Article Number
45	15 – 35 mm	CA4515
57	15 – 40 mm	CA5715



STOMOCUR® Clic Alginat konvex

Plate size 118 mm (Ø), package contents 3 pieces

	Stoma Opening	Plateau Size	Article Number
45	15 – 35 mm	40 mm	CA4515V2
	15 – 25 mm	30 mm	CA4515V3
57	15 – 43 mm	50 mm	CA5715V1
	15 – 35 mm	40 mm	CA5715V2
	15 – 25 mm	30 mm	CA5715V3



STOMOCUR® Select Clic

Window, filter, soft cover both sides, wpouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CC45
57	57 mm	CC57
70	70 mm	CC70



STOMOCUR® Protect Clic small

Window, filter, soft cover both sides, 'Safe and Easy' closure, pouch size 150 x 190 mm, capacity 400 ml, package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CI45S

Soft and discreet closure



STOMOCUR® Protect Clic

Window, filter, soft cover both sides, 'Safe and Easy' closure, pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces

	Locking Ring Size	Article Number
45	45 mm	CI45
57	57 mm	CI57
70	70 mm	CI70

Soft and discreet closure



STOMOCUR® Drainage Protect Clic transparent

Transparent, filter, soft cover one side, large soft outlet, pouch size 150 x 220 mm, capacity 540 ml, package contents 30 pieces

	Locking Ring Size	Article Number
70	70 mm	CDP70T



Our extensive urostomy range is designed for maximum capacity with minimal bulk. The unique anti-twist pouch shape makes emptying easier. Available with soft bung or tap closure for increased safety and convenience, baseplates include flat and convex skin protection.

In addition to our Combiflange hydrocolloid skin protection, our urostomy portfolio includes Microskin baseplates with a semi-permeable adhesive surface. This ultra-thin, elastic and breathable product easily conforms to body contours and folds, significantly reducing the risk of leakage. Microskin is breathable and water vapor-permeable for healthier skin, providing the user ultimate freedom of movement.

Available:


- Three different skin protection options: Combiflange, Microskin, gelatine-free
- Available in flat and convex
- As one-piece and two-piece systems

Our urostomy pouches have the following features:

- Optional 'Safe' soft bung or 'Twist and Safe' tap closure
- Drip cap for 'Twist and Safe' outlet
- Anti-twist shape
- Easy connection to night time drainage systems
- Integral valve to prevent urine backflow
- 'Window' technology for easy positioning
- Soft non-woven cover - quiet, non-cling, comfortable
- Also available with a transparent pouch back



One-Piece Supply Urostomy



STOMOCUR® Safe

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 60 mm	USH10F
15 mm	USH15F
20 mm	USH20F
25 mm	USH25F
30 mm	USH30F

STOMOCUR® Safe transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Article Number
10 – 60 mm	USH10FT



STOMOCUR® Twist & Safe

Window, soft cover both sides, tap valve, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
10 – 60 mm	USH10H
15 mm	USH15H
20 mm	USH20H
25 mm	USH25H
30 mm	USH30H

STOMOCUR® Twist & Safe transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Article Number
10 – 60 mm	USH10HT



Soft and comfortable

STOMOCUR® Safe konvex

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (∅), pouch size 130 × 195 mm, capacity 420 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1F
15 mm		USH15V1F
20 mm		USH20V1F
25 mm		USH25V1F
30 mm		USH30V1F
35 mm		USH35V1F
10 - 35 mm	40 mm	USH10V2F
15 mm		USH15V2F
20 mm		USH20V2F
25 mm		USH25V2F
30 mm		USH30V2F
10 - 25 mm	30 mm	USH10V3F
15 mm		USH15V3F
20 mm		USH20V3F

STOMOCUR® Safe konvex transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1FT
10 - 35 mm	40 mm	USH10V2FT
10 - 25 mm	30 mm	USH10V3FT



'Twist and Safe' tap

STOMOCUR® Twist & Safe konvex



Window, soft cover both sides, tap valve, adhesive surface 118 mm (∅), pouch size 130 × 195 mm, capacity 420 ml, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1H
15 mm		USH15V1H
20 mm		USH20V1H
25 mm		USH25V1H
30 mm		USH30V1H
35 mm		USH35V1H
10 - 35 mm	40 mm	USH10V2H
15 mm		USH15V2H
20 mm		USH20V2H
25 mm		USH25V2H
30 mm		USH30V2H
10 - 25 mm	30 mm	USH10V3H
15 mm		USH15V3H
20 mm		USH20V3H

STOMOCUR® Twist & Safe konvex transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Plateau Size	Article Number
10 - 43 mm	50 mm	USH10V1HT
10 - 35 mm	40 mm	USH10V2HT
10 - 25 mm	30 mm	USH10V3HT



STOMOCUR® Safe medium konvex

Window, soft cover both sides, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 10 pieces



Soft and comfortable

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1F
15 mm		USH15M1F
20 mm		USH20M1F
25 mm		USH25M1F
30 mm		USH30M1F
35 mm		USH35M1F
10 - 40 mm	45 mm	USH10M2F
15 mm		USH15M2F
20 mm		USH20M2F
25 mm		USH25M2F
10 - 28 mm	35 mm	USH10M3F


STOMOCUR® Safe medium konvex transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1FT
10 - 40 mm	45 mm	USH10M2FT
10 - 28 mm	35 mm	USH10M3FT

STOMOCUR® Twist & Safe medium konvex

Window, soft cover both sides, tap valve, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 10 pieces



'Twist and Safe' tap

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1H
15 mm		USH15M1H
20 mm		USH20M1H
25 mm		USH25M1H
30 mm		USH30M1H
35 mm		USH35M1H
10 - 40 mm	45 mm	USH10M2H
15 mm		USH15M2H
20 mm		USH20M2H
25 mm		USH25M2H
10 - 28 mm	35 mm	USH10M3H

STOMOCUR® Twist & Safe medium konvex transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Plateau Size	Article Number
10 - 47 mm	55 mm	USH10M1HT
10 - 40 mm	45 mm	USH10M2HT
10 - 28 mm	35 mm	USH10M3HT



STOMOCUR® Safe Clic

Window, soft cover both sides, soft outlet, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Locking Ring Size	Article Number
45 mm	CU45F
57 mm	CU57F

STOMOCUR® Safe Clic transparent

Transparent, soft cover one side, soft outlet

Locking Ring Size	Article Number
45 mm	CU45FT
57 mm	CU57FT



STOMOCUR® Safe Microskin

Window, soft cover both sides, soft outlet, adhesive surface 135 x 100 mm, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
13 - 50 mm	USM13F

STOMOCUR® Safe Microskin transparent

Transparent, soft cover one side, soft outlet

Stoma Opening	Article Number
13 - 50 mm	USM13FT



STOMOCUR® Twist & Safe Clic

Window, soft cover both sides, tap valve, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Locking Ring Size	Article Number
45 mm	CU45H
57 mm	CU57H

STOMOCUR® Twist & Safe Clic transparent

Transparent, soft cover one side, tap valve

Locking Ring Size	Article Number
45 mm	CU45HT
57 mm	CU57HT



STOMOCUR® Twist & Safe Microskin

Window, soft cover both sides, tap valve, adhesive surface 135 x 100 mm, pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
13 - 50 mm	USM13H

STOMOCUR® Twist & Safe Microskin transparent

Transparent, soft cover one side, tap valve

Stoma Opening	Article Number
13 - 50 mm	USM13HT



To select the suitable baseplates, please see page 24 (STOMOCUR® Clic) and page 25 (STOMOCUR® Clic konvex).

One-Piece Supply Urostomy



STOMOCUR® Safe gelatinefrei transparent

Transparent, soft cover one side, soft outlet, adhesive surface 118 mm (Ø), pouch size 130 x 195 mm, capacity 420 ml, package contents 30 pieces

Stoma Opening	Article Number
15 – 60 mm	USG15FT



Soft and comfortable



Accessories



STOMOCUR® Uro Adapter

For connecting bed bags to the STOMOCUR® urostomy pouches, package contents 5 pieces

Article Number
USZ1



STOMOCUR® Drip Cap

For use with the STOMOCUR® 'Twist and Safe' tap, package contents 5 pieces

Article Number
USZ2



gelatinefrei

For ostomy patients who do not use products with gelatine or other animal-based ingredients, we offer our pouch assortment with gelatine-free skin protection. By using natural, plant-based ingredients, these pouches are equally skin-friendly and do not compromise on comfort and safety.

The skin protection used is made of hypoallergenic hydrocolloid without gelatine and is available in our one-piece closed Select colostomy pouches, one-piece drainable Protect ileostomy pouches and Safe urostomy pouches.

Features:

- Thin and flexible
- Particularly gentle on the skin
- Strong adhesive



STOMOCUR® Select gelatinefrei

Window, filter, soft cover both sides, adhesive surface 110 mm (Ø), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	CSG20



Soft and discreet closure

STOMOCUR® Protect gelatinefrei

Window, filter, soft cover both sides, 'Safe and Easy' closure, adhesive surface 110 mm (Ø), pouch size 150 x 220 mm, capacity 530 ml, package contents 30 pieces

Stoma Opening	Article Number
20 – 60 mm	IPG20



SPECIAL SUPPLY

Our comprehensive range of products for specialised care has been developed for medical professionals and includes various types of drainage pouches as well as pouches for paediatric care.

Paediatric Care:

- Two options for premature babies up to 900 g and infants up to 5000 g
- Gentle and skin-friendly hydrocolloid
- Soft non-woven cover on one side, transparent back for easy observation
- Closed pouch can be cut open and used as a drainable pouch with a clip or an adhesive tape

Drainage Protect Pouch:

- Available in flat, convex and medium convex options
- Soft non-woven cover with a window for easy positioning and observation
- MPC filter for odour control
- High output drain with bung closure
- Hypoallergenic hydrocolloid as skin protection

Drainage Pouches:

- Large, flexible skin protection plates combined with flexible, transparent pouches – ideal for promoting healing with easy observation
- High-output drainage ports with soft outlet
- Integrated access openings allow easy access to the wound without removing the pouch
- A range of shapes and sizes to meet the patients' needs

STOMOCUR® Safe gelatinefrei transparent - Flat urostomy pouches with gelatine-free skin protection can be found on page 36.



For premature babies up to 900 g

STOMOCUR® Premature

Transparent, filter, soft cover one side, adhesive surface 40 mm, pouch size 110 mm x 45 mm, capacity 20 ml, package contents 20 pieces

Stoma Opening	Article Number
3 mm	P3003



For infants up to 5000 g

STOMOCUR® Premature Combi

Transparent, soft cover one side, adhesive surface 40 mm, pouch size 210 mm x 45 mm, capacity 90 ml, package contents 30 pieces

Stoma Opening	Article Number
3 mm	P3203



STOMOCUR® Drainage Extra 'Access Port'

Extra Small, opening 10 mm, foil: transparent, package contents 1 piece

Article Number	Max. Cutability	Access Opening
SDE2-XSA	90 x 180 mm	ø 120 mm



STOMOCUR® Drainage Extra 'Access Port'

Small, opening 10 mm, foil: transparent, package contents 1 piece

Article Number	Max. Cutability	Access Opening
SDE2-SA	240 x 155 mm	ø 120 mm



STOMOCUR® Drainage Extra 'Access Port'

Medium, opening 10 mm, foil: transparent, package contents 1 piece

Article Number	Max. Cutability	Access Opening
SDE2-MA	85 x 250 mm	ø 120 mm



STOMOCUR® Drainage Extra 'Access Port'

Large, opening 10 mm, foil: transparent, package contents 1 piece

Article Number	Max. Cutability	Access Opening
SDE2-LA	160 x 200 mm	ø 120 mm

With 2 drainage outlets



STOMOCUR® Drainage Extra 'Access Port'

Extra Large, opening 10 mm, foil: transparent, package contents 1 piece

Article Number	Max. Cutability	Access opening
SDE2-XXLA	230 x 380 mm	ø 120 mm

Soft and flexible



STOMOCUR® Drainage Protect

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, package contents 10 pieces

Stoma Opening	Article Number
20 - 60 mm	DP5020

STOMOCUR® Drainage small - large

Transparent with printed scale, backstop, exhaust valve with adapter, package contents 30 pieces

Stoma Opening	Capacity	Article Number
3 mm	50 ml	DP2503XS
5 mm	150 ml	DP4005S
5 mm	300 ml	DP4005



STOMOCUR® Drainage Protect konvex

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, convex height 7 mm, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 43 mm	50 mm	DP20V1
15 - 35 mm	40 mm	DP15V2
10 - 25 mm	30 mm	DP10V3

Available in 3 plateau sizes



STOMOCUR® Drainage Protect medium konvex

Window, filter, soft cover both sides, large soft outlet, adhesive surface 118 mm (Ø), pouch size 150 x 220 mm, capacity 540 ml, convex height 3,5 mm, package contents 10 pieces

Stoma Opening	Plateau Size	Article Number
20 - 47 mm	55 mm	DP20M1
15 - 40 mm	45 mm	DP15M2
10 - 28 mm	35 mm	DP10M3

Available in 3 plateau sizes




The comprehensive FORLIFE range of ostomy care products and accessories provides additional safety and reassurance.

- Hydrocolloid skin protection plates
- Skin protection rings
- Nourishing lotion
- Skin protection paste
- Skin protection powder
- Adhesive remover spray
- Superabsorbent gelling agents
- Fixing strips
- Ostomy belts
- Urostomy adapters for connecting bed and leg bags
- Urostomy drip cap
- Scissors



Accessories



STOMOCUR® Skin Protection Plates

With self-adhesive hydrocolloid skin protection compound, material thickness 0,6 mm, package contents 5 pieces

Size	Article Number
100 x 100 mm	HP101006
150 x 150 mm	HP151506
200 x 200 mm	HP202006
370 x 220 mm	HP372206



STOMOCUR® Skin Protection Plates mouldable

With self-adhesive hydrocolloid skin protection compound, can be individually cut and shaped, material thickness 2,4 mm, package contents 1 piece

Size	Article Number
220 x 220 mm	HP222224



STOMOCUR® Skin Protection Rings

With hydrocolloid skin protection material to seal and to protect the skin around the stoma, outer diameter 48 mm or 63 mm, material thickness 2,4 mm, package contents 30 pieces

Size	Article Number
10 mm	48101
15 mm	48151
20 mm	48201
25 mm	63251
30 mm	63301
35 mm	63351
40 mm	63401



Accessories



STOMOCUR® Fixing Strips

Skin-friendly and highly elastic, extends the wearing time, prevents the baseplate from detaching and rolling up, more safety and comfort, package contents 20 pieces

Article Number
FS20



STOMOCUR® Skin Protection Paste

For external use on the skin around the stoma, fills unevenness between the stoma and the skin protection plate, available with and without alcohol, package contents 1 piece

Article Number	Contents
CP2050 (with alcohol)	56 g
PASTE0A (without alcohol)	56 g



STOMOCUR® Skin Protection Powder

The powder absorbs the moisture of the skin around the stoma, ensuring better adhesion of the stoma supply, package contents 1 piece

Article Number	Contents
HP30	30 g



Accessories



STOMOCUR® Safety

The skin care lotion is applied very thinly to the skin around the stoma to prevent skin irritation, package contents 1 piece

Article Number	Contents
SSAFETY100	100 ml



STOMOCUR® EmplaSectal

For easy and skin-friendly removal of the skin protection plates, handy as a spray, package contents 1 piece

Article Number	Contents
EMPLASECTAL	50 ml



STOMOCUR® Ostomy Belt

Storage in two lengths, adjustable, width 3 cm, package contents 1 piece

Article Number	Length
AB01	100 cm
AB02	120 cm



STOMOCUR® SGX

Fast-acting, super-absorbent gelling agent, easy to use, helpful for ileostomists and colostomists with recipients with more fluid output, reduces pouch swelling, pouch noise and odour, ideal for travel, package contents 30 pieces



Article Number
SGX30



Accessories



STOMOCUR® Uro Adapter

For connecting bed bags to the STOMOCUR® urostomy pouches, package contents 5 pieces

Article Number
USZ1



STOMOCUR® Uro Drip Cap

To attach the STOMOCUR® 'Twist and Safe' valve, package contents 5 pieces

Article Number
USZ2



STOMOCUR® Scissors

Scissors with rounded blades for cutting the baseplate openings, package contents 1 piece

Article Number
SSchere

Table of contents by item number

Article Number	Detail	Page
48101	10 mm	45
48151	15 mm	45
48201	20 mm	45
63251	25 mm	45
63301	30 mm	45
63351	35 mm	45
63401	40 mm	45
AB01	100 cm	47
AB02	120 cm	47
CA4515V2	15 mm - 35 mm	25
CA4515V3	15 mm - 25 mm	25
CA5715V1	15 mm - 43 mm	25
CA5715V2	15 mm - 35 mm	25
CA5715V3	15 mm - 25 mm	25
CA4515	15 mm - 35 mm	24
CA5715	15 mm - 40 mm	24
CC45	45 mm	26
CC57	57 mm	26
CC70	70 mm	26
CDP70T	70 mm	27
CI45	45 mm	26
CI45S	45 mm	27
CI57	57 mm	26
CI70	70 mm	26
CK4515	15 mm - 35 mm	24
CK4515V2	15 mm - 35 mm	25
CK4515V3	15 mm - 25 mm	25
CK4520	20 mm	24
CK4520V2	20 mm	25
CK4525	25 mm	24
CK4525V2	25 mm	25
CK4530	30 mm	24
CK4530V2	30 mm	25
CK5715	15 - 40 mm	24
CK5715V1	15 - 43 mm	25
CK5715V2	15 - 35 mm	25
CK5715V3	15 - 25 mm	25
CK5720V1	20 mm	25
CK5725V1	25 mm	25
CK5730V1	30 mm	25
CK5735	35 mm	24
CK5735V1	35 mm	25
CK5740	40 mm	24
CK7015	15 mm - 50 mm	24
CK7045	45 mm	24
CK7050	50 mm	24
CP2050	56 g	46
CSA10V3	10 mm - 25 mm	15
CSA15V2	15 mm - 35 mm	15
CSA15V3	15 mm	15
CSA20	20 mm - 60 mm	21
CSA20S	20 mm - 60 mm	21
CSA20V1	20 mm - 43 mm	15
CSA20V2	20 mm	15
CSA20V3	20 mm	15
CSA25	25 mm	21

Table of contents by item number

Article Number	Detail	Page
CSA25S	25 mm	21
CSA25V1	25 mm	15
CSA25V2	25 mm	15
CSA30	30 mm	21
CSA30S	30 mm	21
CSA30V1	30 mm	15
CSA30V2	30 mm	15
CSA35	35 mm	21
CSA35S	35 mm	21
CSA35V1	35 mm	15
CSA40	40 mm	21
CSA40S	40 mm	21
CSA40V1	40 mm	15
CSC20	20 mm – 50 mm	10
CSG20	20 mm – 60 mm	38
CSH10M3	10 mm – 28 mm	18
CSH10V3	10 mm – 25 mm	14
CSH10XS	10 mm – 50 mm	10
CSH10XST	10 mm – 50 mm	10
CSH15M2	15 mm – 40 mm	18
CSH15V2	15 mm – 35 mm	14
CSH15V3	15 mm	14
CSH20	20 mm – 60 mm	9
CSH20L	20 mm – 90 mm	9
CSH20M1	20 mm – 47 mm	18
CSH20M2	20 mm	18
CSH20S	20 mm – 60 mm	10
CSH20V0	20 mm – 57 mm	14
CSH20V1	20 mm – 43 mm	14
CSH20V2	20 mm	14
CSH20V3	20 mm	14
CSH25	25 mm	9
CSH25L	25 mm	9
CSH25M1	25 mm	18
CSH25M2	25 mm	18
CSH25S	25 mm	10
CSH25V1	25 mm	14
CSH25V2	25 mm	14
CSH30	30 mm	9
CSH30L	30 mm	9
CSH30M1	30 mm	18
CSH30S	30 mm	10
CSH30V1	30 mm	14
CSH30V2	30 mm	14
CSH35	35 mm	9
CSH35L	35 mm	9
CSH35M1	35 mm	18
CSH35S	35 mm	10
CSH35V1	35 mm	14
CSH40	40 mm	9
CSH40L	40 mm	9
CSH40M1	40 mm	18
CSH40S	40 mm	10
CSH40V1	40 mm	14
CU45F	45 mm	34
CU45FT	45 mm	34

Article Number	Detail	Page
CU45H	45 mm	34
CU45HT	45 mm	34
CU57F	57 mm	34
CU57FT	57 mm	34
CU57H	57 mm	34
CU57HT	57 mm	34
DP10M3	10 mm – 28 mm	19 u. 43
DP10V3	10 mm – 25 mm	16 u. 43
DP15M2	15 mm – 40 mm	19 u. 43
DP15V2	15 mm – 35 mm	16 u. 43
DP20M1	20 mm – 47 mm	19 u. 43
DP20V1	20 mm – 43 mm	16 u. 43
DP2503XS	3 mm	42
DP4005	5 mm	42
DP4005S	5 mm	42
DP5020	20 mm – 60 mm	12 u. 43
EMPLASECTAL	50 ml	47
FS20		46
HP101006	100 × 100 mm	45
HP151506	150 × 150 mm	45
HP202006	200 × 200 mm	45
HP222224	220 × 220 mm	45
HP30	30 g	46
HP372206	370 × 220 mm	45
IPH15T	15 mm – 60 mm	11
IPA10V3	10 mm – 25 mm	15
IPA15	15 mm – 60 mm	22
IPA15V1	15 mm – 43 mm	15
IPA15V2	15 mm – 35 mm	15
IPA20	20 mm	22
IPA20V1	20 mm	15
IPA25	25 mm	22
IPA25V1	25 mm	15
IPA30	30 mm	22
IPA30V1	30 mm	15
IPA35	35 mm	22
IPA35V1	35 mm	15
IPA40	40 mm	22
IPA40V1	40 mm	15
IPG20	20 mm – 60 mm	38
IPH10M3	10 mm – 28 mm	18
IPH10M3T	10 mm – 28 mm	18
IPH10V3	10 mm – 25 mm	14
IPH15M1	15 mm – 47 mm	18
IPH15M1T	15 mm – 47 mm	18
IPH15M2	15 mm – 40 mm	18
IPH15M2T	15 mm – 40 mm	18
IPH15T	15 mm – 60 mm	11
IPH15S	15 mm – 60 mm	11
IPH15V1	15 mm – 43 mm	14
IPH15V2	15 mm – 35 mm	14
IPH15V3	15 mm	14
IPH20	20 mm	11
IPH20M1	20 mm	18
IPH20M2	20 mm	18
IPH20V0	20 mm – 57 mm	14

Table of contents by item number

Article Number	Detail	Page
IPH20V1	20 mm	14
IPH20V2	20 mm	14
IPH20V3	20 mm	14
IPH2070	20 mm - 70 mm	11
IPH25	25 mm	11
IPH25M1	25 mm	18
IPH25M2	25 mm	18
IPH25V1	25 mm	14
IPH25V2	25 mm	14
IPH30	30 mm	11
IPH30M1	30 mm	18
IPH30V1	30 mm	14
IPH30V2	30 mm	14
IPH35	35 mm	11
IPH35M1	35 mm	18
IPH35V1	35 mm	14
IPH40	40 mm	11
IPH40M1	40 mm	18
IPH40V1	40 mm	14
P3003	3 mm	40
P3203	3 mm	40
PASTEOA	56 g	46
SDE2-LA	160 x 200 mm	41
SDE2-MA	85 x 250 mm	41
SDE2-SA	240 x 155 mm	41
SDE2-XSA	90 x 180 mm	41
SDE2-XXLA	230 x 380 mm	42
SGX30		47
SSAFETY100	100 ml	47
SSchere		48
USG15FT	15 mm - 60 mm	36
USH10F	10 mm - 60 mm	29
USH10FT	10 mm - 60 mm	29
USH10H	15 mm - 60 mm	29
USH10HT	10 mm - 60 mm	29
USH10M1F	10 mm - 47 mm	32
USH10M1FT	10 mm - 47 mm	32
USH10M1H	10 mm - 47 mm	33
USH10M1HT	10 mm - 47 mm	33
USH10M2F	10 mm - 40 mm	32
USH10M2FT	10 mm - 40 mm	32
USH10M2H	10 mm - 40 mm	33
USH10M2HT	10 mm - 40 mm	33
USH10M3F	10 mm - 28 mm	32
USH10M3FT	10 mm - 28 mm	32
USH10M3H	10 mm - 28 mm	33
USH10M3HT	10 mm - 28 mm	33
USH10V1F	10 mm - 43 mm	30
USH10V1FT	10 mm - 43 mm	30
USH10V1H	10 mm - 43 mm	31
USH10V1HT	10 mm - 43 mm	31
USH10V2F	10 mm - 35 mm	30
USH10V2FT	10 mm - 35 mm	30
USH10V2H	10 mm - 35 mm	31
USH10V2HT	10 mm - 35 mm	31
USH10V3F	10 mm - 25 mm	30

Article Number	Detail	Page
USH10V3FT	10 mm - 25 mm	30
USH10V3H	10 mm - 25 mm	31
USH10V3HT	10 mm - 25 mm	31
USH15F	15 mm	29
USH15H	15 mm	29
USH15M1F	15 mm	32
USH15M1H	15 mm	33
USH15M2F	15 mm	32
USH15M2H	15 mm	33
USH15V1F	15 mm	30
USH15V1H	15 mm	31
USH15V2F	15 mm	30
USH15V2H	15 mm	31
USH15V3F	15 mm	30
USH15V3H	15 mm	31
USH20F	20 mm	29
USH20H	20 mm	29
USH20M1F	20 mm	32
USH20M1H	20 mm	33
USH20M2F	20 mm	32
USH20M2H	20 mm	33
USH20V1F	20 mm	30
USH20V1H	20 mm	31
USH20V2F	20 mm	30
USH20V2H	20 mm	31
USH20V3F	20 mm	30
USH20V3H	20 mm	31
USH25F	25 mm	29
USH25H	25 mm	29
USH25M1F	25 mm	32
USH25M1H	25 mm	33
USH25M2F	25 mm	32
USH25M2H	25 mm	33
USH25V1F	25 mm	30
USH25V1H	25 mm	31
USH25V2F	25 mm	30
USH25V2H	25 mm	31
USH30F	30 mm	29
USH30H	30 mm	29
USH30M1F	30 mm	32
USH30M1H	30 mm	33
USH30V1F	30 mm	30
USH30V1H	30 mm	31
USH30V2F	30 mm	30
USH30V2H	30 mm	31
USH35M1F	35 mm	32
USH35M1H	35 mm	33
USH35V1F	35 mm	30
USH35V1H	35 mm	31
USM13F	13 mm - 50 mm	35
USM13FT	13 mm - 50 mm	35
USM13H	13 mm - 50 mm	35
USM13HT	13 mm - 50 mm	35
USZ1		36 and 48
USZ2		36 and 48


Your contact details:





For Life
Produktions- und Vertriebsgesellschaft
für Heil- und Hilfsmittel mbH

Am Studio 16
12489 Berlin, Germany


mail@forlife.info
www.forlife.info

 @For Life Produktions- und Vertriebsgesellschaft für Heil- und Hilfsmittel mbH

 @For Life Produktions- und Vertriebsgesellschaft für Heil- und Hilfsmittel mbH

 @For Life Berlin

 @FORLIFE_BERLIN

 @forlife_deutschland


BERLIN

