


KONTICUR®

Urostomy Bed Bag KUBB2000



Safety and Comfort for a Restful Sleep at Night

The KONTICUR® Urostomy Bed Bag can be connected to all FORLIFE urostomy pouches easily thanks to the integrated adapter. The kink-resistant tube enables smooth drainage and the one-hand drain can be operated safely even with limited motor skills.

- » Integrated adapter can be connected to all FORLIFE urostomy pouches
- » Easy to operate tap outlet designed for one-hand operation
- » 120 cm long kink-resistant tube

Order information:

Item name	Item number	Tube	Capacity	Package contents
KONTICUR® Urostomy Bed Bag	KUBB2000	120 cm	2000 ml	5 pieces

